

# URINSAFE™

Reagent Strips for Urinalysis

pH	Meaning
5.0 to 7.0	Normal
7.0 to 8.0	Normal - add .005 to SG result
8.0	<b>Caution*</b>

Specific Gravity	Hydration Status
1.000 to 1.005	<b>Excellent **</b>
1.005 to 1.010	Very Good
1.010 to 1.015	Good
1.015 to 1.020	Fair
1.020 to 1.025	Caution - Borderline
>1.025	Poor - Unsatisfactory

**100 Strips**




Call us on 1800 681 544 or visit [www.lfafirstresponse.com.au](http://www.lfafirstresponse.com.au)

## TESTING AND READING TIME

Rev. 04/2019

<b>pH 60s</b>	5.0	6.0	6.5	7.0	7.5	8.0	8.5
<b>Specific Gravity 45s</b>	1.000	1.005	1.010	1.015	1.020	1.025	1.030



\* A pH of 8.0 **may** indicate a urinary tract infection or an old sample. Retest with a fresh sample. If the pH is still 8.0 seek medical advice.

\*\* With persistent low readings, **approaching an SG of 1.000**, seek medical advice as it is necessary to rule out over hydration or dilutional hyponatremia.



Distributed by  
**LFA First Response**  
 Australian helpline 1800 681 544  
[www.lfafirstresponse.com.au](http://www.lfafirstresponse.com.au)